

# Dynamic RF Muscle Toning System - Official Clinical Overview & Technical Datasheet

## DEVICE IDENTIFICATION AND CLINICAL PURPOSE

The Dynamic RF Muscle Toning System is a class IIa medical aesthetic device indicated for non-invasive muscle stimulation, augmentation of muscle tone, and improvement of abdominal and gluteal definition. The system delivers simultaneous multi-channel radiofrequency energy combined with synchronized electrical muscle stimulation (EMS) to induce supramaximal muscle contractions while promoting thermal collagen modulation in overlying dermal tissues.



## INTERNAL HARDWARE TOPOLOGY

The system architecture integrates a digitally controlled RF generator (1.0 MHz carrier frequency) with a dedicated EMS waveform synthesizer operating independently across four output channels. Each channel is managed by a galvanically isolated power stage utilizing silicon carbide MOSFET switching technology, ensuring sub-millisecond response to impedance fluctuations. The master control unit employs a dual-core real-time processor implementing adaptive feedback algorithms that monitor tissue impedance (range: 20–600 ohms) and dynamically modulate RF power delivery to maintain preset energy targets. A redundant safety interlock system passively disables all outputs when contact integrity falls below specified thresholds.

#### EPIDERMAL PROTECTION MECHANISMS

The system incorporates a three-tier thermal management protocol. First, integrated thermistors embedded within each treatment pad provide real-time surface temperature monitoring with a sampling rate of 50 Hz. Second, a closed-loop circulating water cooling circuit maintains the pad baseplate at  $15^{\circ}\text{C} \pm 2^{\circ}\text{C}$ , dissipating heat generated by RF energy and preventing unintended thermal injury. Third, software-enforced treatment limits cap cumulative RF fluence at 45 kJ per anatomical zone per session. Automatic output termination activates if skin temperature exceeds  $42^{\circ}\text{C}$  or if impedance deviates beyond safe operating windows.

## TREATMENT ADVANTAGES

Unlike conventional EMS devices that elicit only superficial motor unit recruitment, Dynamic RF Muscle Toning achieves volumetric muscle fiber activation up to 100% of the targeted motor unit pool. Concurrent RF energy preconditioning reduces muscle viscosity, enabling deeper and more efficient electrical field penetration. Clinical studies demonstrate an average 16% increase in rectus abdominis muscle thickness after eight biweekly treatments, with corresponding reductions in subcutaneous adipose thickness (12% mean) measured via ultrasound. Patient-reported outcomes indicate an 89% satisfaction rate concerning visible abdominal contour improvement at 12-week post-treatment follow-up.

## SPECIFICATION MATRIX

<b>Parameter</b>	<b>Specification</b>
RF Carrier Frequency	1.0 MHz $\pm$ 5%
EMS Frequency Range	10 Hz – 100 Hz (programmable)
Maximum RF Power per Channel	120 W (continuous), 240 W (peak)
Maximum EMS Current per Channel	120 mA (peak-to-peak)
Number of Independent Channels	4 (each supports RF + EMS)

	simultaneous or sequential)
Impedance Monitoring Range	20 – 600 ohms, adaptive power adjustment
Cooling Mechanism	Thermoelectric + closed-loop water circulation (15 ° C ± 2 ° C pad temperature)
Temperature Safety Limits	Automatic cutoff at 42°C skin surface
Treatment Pad Material	Medical-grade silicone with hydrogel conductive layer
Display Interface	10.1-inch capacitive touchscreen, 1280 x 800 resolution
Power Supply	100–240 VAC, 50/60 Hz, 700 VA
Dimensions (Main Unit)	450 mm (W) x 400 mm (D) x 1100 mm (H)
Weight	28.5 kg (main console)
Standards Compliance	CE MDR (EU 2017/745), FDA 510(k) K221345, ISO 13485, IEC 60601 series

## REGULATORY COMPLIANCE

The Dynamic RF Muscle Toning System carries CE marking under Medical

Device Regulation (EU) 2017/745 and holds FDA 510(k) clearance (K221345) for muscle toning and body contouring adjunctive indications. Manufacturing facilities are ISO 13485:2016 certified. Electrical safety compliance includes IEC 60601-1:2012 (medical electrical equipment), IEC 60601-2-2:2017 (high-frequency surgical equipment and RF therapy devices), and IEC 60601-1-11:2015 (home healthcare environment requirements). EMC conformity to IEC 60601-1-2:2014 Edition 4.1 is verified. Biocompatibility of all patient-contacting components (pads, straps, gel interfaces) meets ISO 10993-5 (cytotoxicity) and ISO 10993-10 (irritation and skin sensitization) standards.

