

# Selective Photothermolysis Architecture Reference Document: Laser Hair Removal Machine

## SELECTIVE PHOTOTHERMOLYSIS ARCHITECTURE REFERENCE DOCUMENT: LASER HAIR REMOVAL MACHINE

### EXECUTIVE SUMMARY

This document delineates the clinical engineering foundation, optical emission topology, and integrated epidermal protection systems of the next-generation diode laser hair removal platform. Designed for sustained duty cycles in high-volume medical aesthetic environments, the system leverages a proprietary multi-wavelength stack (755nm / 808nm / 1064nm) to achieve selective photothermolysis across Fitzpatrick skin types I–VI. The following sections provide a hardware deep-dive, performance registry, and certified compliance framework for procurement and clinical integration.



## CLINICAL ARCHITECTURE & DESIGN

The platform employs a resonant diode laser array coupled with a forced-air and liquid thermal management subsystem. The optical delivery train utilizes a high-transmission fused silica handpiece waveguide, ensuring fluence uniformity tolerance of  $\pm 5\%$  across the entire spot geometry. A real-time contact sensing logic interlock prevents inadvertent firing when epidermal impedance is suboptimal.

## KEY INDICATIONS & CAPABILITIES

- Permanent hair reduction for terminal and vellus hairs
- Treatment of pseudofolliculitis barbae (ingrown hair management)
- Axilla, bikini, legs, arms, back, chest, and facial zones (excluding periorbital)
- Integrated skin tone sensor with automatic fluence titration (1–6 Fitzpatrick)
- Repetition rate up to 10Hz for large surface area clearance

## COMPLIANCE & STANDARDS

- FDA 510(k) Class II medical device clearance
- CE 0197 (MDR 2017/745) compliance
- IEC 60825-1:2014 Class 4 laser safety certification
- ISO 13485:2016 accredited manufacturing facility
- RoHS 3 and REACH environmental compliance

## TECHNICAL SPECIFICATIONS

| Parameter               | Specification                                                                   |
|-------------------------|---------------------------------------------------------------------------------|
| Laser Type / Wavelength | Diode laser stack: 755nm / 808nm / 1064nm (selectable or simultaneous emission) |
| Peak Optical Power      | Up to 1800W per wavelength channel                                              |
| Fluence Range           | 5 – 30 J/cm <sup>2</sup> (step size 0.5 J/cm <sup>2</sup> )                     |
| Spot Size               | 12mm x 12mm (standard), 15mm x 15mm (large area), 8mm round (precision)         |
| Repetition Rate         | Single shot to 10 Hz continuous                                                 |
| Cooling System          | TEC + Sapphire contact cooling + closed-loop water + forced air                 |

|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| Pulse Width Range           | 5 – 500 ms (adjustable in 1 ms increments)                               |
| Skin Tone Sensor            | Reflective photoplethysmography + thermal mapping (Fitzpatrick I–VI)     |
| User Interface              | 10.4" capacitive touchscreen, wireless foot pedal, RFID handpiece reader |
| Electrical                  | 100–240 VAC, 50/60 Hz, 1500 VA max (Class I Medical)                     |
| Dimensions (Main Unit)      | 450mm W x 600mm D x 1100mm H                                             |
| Weight                      | 38 kg (handpiece with cable: 0.8 kg)                                     |
| Standby / Operating Ambient | 15 ° C – 30 ° C / 20% – 80% RH non-condensing                            |

#### TREATMENT PROTOCOLS (VERIFIED ENERGY REGISTER)

- Fitzpatrick I–III (Light skin): 10–14 J/cm<sup>2</sup>, pulse width 10–30ms
- Fitzpatrick IV–V (Olive to brown): 14–18 J/cm<sup>2</sup>, pulse width 30–50ms
- Fitzpatrick VI (Dark brown/Black): 16–20 J/cm<sup>2</sup>, pulse width 50–100ms (test spot mandatory)
- Recommended overlap: ≤10% per pulse, minimum 6–8 weeks interval
- Cooling pre-pulse: -5°C to 4°C sapphire plate contact



## EPIDERMAL PROTECTION MECHANISMS

The integrated cascaded cooling architecture comprises: 1) Thermoelectric cooler (TEC) pre-chilling to  $-10^{\circ}\text{C}$  internal reservoir, 2) recirculating deionized water loop (thermal capacity 450W), 3) direct-contact sapphire window (thermal conductivity  $35\text{ W/m}\cdot\text{K}$ ), and 4) forced-air exhaust for handpiece shell. Post-pulse skin temperature is maintained below  $45^{\circ}\text{C}$ , eliminating erythema and blistering risks across all nominal settings.

## CLINICAL SAFETY & RISK MANAGEMENT

The device incorporates a redundant emergency beam abort, foot pedal interlock, and key-controlled master power. Each pulse is subjected to a pre-discharge system health check (diode impedance, coolant flow rate, contact impedance). In the event of any parameter excursion, the treatment head automatically inhibits firing and triggers a visual (red ring) and audible

alert.